

# Oxivir Excel® CE Wipe



Bactericidal



Virucidal



Yeastcidal



- Ready to use for faster disinfection
- Disposable wipe for ultra hygiene
- Convenient, no preparation or laundry required
- Kills MDRO's such as MRSA, CRE, VRE, CRPSA
- Effective against Bacteria, Viruses and Yeast in as little as 30 seconds
- Full virucidal claim including non-enveloped viruses such as norovirus & poliovirus
- Extensively tested against the new European norms
- Carries yeastcidal claim against Candida albicans
- Powerful surfactant package delivers excellent cleaning performance including descaling and lime limitation
- Excellent user safety profile reduces worker risk and absence
- No hazard warnings in dilution
- Does not contain respiratory irritants
- Does not contain salicylic acid
- AHP formulation is gentle on surfaces extending equipment life
- CEIIa certification allows use on for non-invasive medical devices
- Available as a concentrated, RTU or wipe
- Suitable for use in high risk areas including ICU, OR, NICU and patient rooms
- Applied for VAH listing

## Surface Compatibility



Oxivir Excel Wipes can be used on any hard surface usually encountered in a healthcare environment including stainless steel, glass, acrylics and PVC. It is not suitable for use on marble, copper, brass or other surfaces sensitive to oxidative agents.

Efficacy	Test Method	Contact Time Undiluted	Dirty/Clean Conditions
Virucidal (BVD & Mod VVA)	DVV limited virucidal test	30 sec	Dirty (FCS)
Virucidal (Adeno, MNV, Polio)	EN 14476 fully virucidal statement	30 sec	Dirty
Bactericidal	EN 1276	1 min	Dirty
Yeastcidal	EN 13624	5 min	Dirty
Bactericidal	EN 13727	30 sec	Dirty
Yeastcidal	EN 1650	5 min	Dirty
Bactericidal Yeastcidal	EN 16615	1 min	Dirty

Test reports are available on request

## Authorized Representative & Distributor:

**K-One MediTech Sdn Bhd**

68, Lengkok Rishah 2, Kawasan Perindustrian Silibin, 30100 Ipoh, Perak, Malaysia

[www.k-one.com](http://www.k-one.com)